

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

IL-17RD Fusion Protein



Basic Information

Catalog Number:

Ag17286

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

ALDLWEDFSLCREGQREWWIQKIHESQFIIVVCSKGM
KYFVDKKNYKHKGGGRGSGKGLFLVAVSAIAEKLKQ
AKQSSSAALSKFIAVYFDYSCEGDVPGILDSTKYRLM
DNLPLQCSHLHSRDHGLQEPGQHTROGSRNNYFRSK
SGRSLYVAICNMHQFIDEEDWFEKQFVPHPPPLRY
REPVLEKFDSGLVLDVDMCKPGPESDFCLKVEAAVLG
ATGPADSQHESQHGGLDQDGEARPALDGSAAALQPLL
HTVKA GSPDMPRDSGIYDSSVPSSELSPLMEGLSTD
QTETSSLTESVSSSSGLGEEPPALPSKLLSSGCKADL
GCRSYDELHAVAPL

(358-707 aa encoded by BC111702)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.