FOR IN VITRO RESEARCH USE ONLY. NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human TGFA protein



Basic Information

Catalog Number:

Ag17312

Size:

50 μg Form:

Available lyophilized

Species: human

Expression Source:

 $\it e\,coli.$ -derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

ENSTSPLSDPPVAAAVVSHFNDCPDSHTQFCFHGTC RFLVQEDKPACVCHSGYVGARCEHADLLAVVAASQ KKQAITALV

(24-103 aa encoded by BC005308)

Selected Validation Data