

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human MACC1 protein



## Basic Information

Catalog Number:

Ag17445

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

ITTPDPTPNLKRSLNPGYLQKKEEIKSAPLSPKILVKYP  
TFQDKTLNFSNYGVTLKAVLRQSKI DYFLEYFKGDTIA  
LLGEGKVKAIGQSKVKEWYV GVLRGKIGLVHCKNVKV  
ISKEQVMFMSDSVFTTRNLEQIVLPLKLTYYISVVLV  
LVSEKVYDWKVLADVLGYSHLSLEDFDQIQADKESEK  
VSYVIKKLKEDCHTERNTRKFLYELIVALLKMDQCQELV  
ARLIQEAAVLTS AVKLGKGWRELAEKLVRLTKQQMEA  
YEIPHRGNTGDVA VEMMWKPAYDFLYTWSAHYGNN  
YRDVLQDLQSA LDRMKNPVTKHWRELTGVLLVNSLE  
VLRVTAFTSEEV

(502-852 aa encoded by BC137090)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.