

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human DLL1- Specific protein



Basic Information

Catalog Number:

Ag17553

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MTCADGPCFNGGRCSDSPDGGYSCRCVPVGYSGFNC
EKKIDYCSSSPCSNGAKCVDLGDAYLCRCQAGF5GR
HCDDNVDDCASSPCANGGTCRDGVNDF5CTCPPGY
TGRNCSAPVSRCEHAPCHNGATCHERGHRYVCECA
RGYGGPNCQFLLPELPPGPAVVDL5TEKLEGGGPF5PW
VAVCAGVILVLM5LLGCAAVV5CVRLRLQK5HRPPADP
CRGETETM5NLANCQ5REKDISV5IIGATQ5IKNTN5KKA
DFHGDH5ADKNGFK5ARYPAVDYN5LVQDLK5GDDTAV
RDAHSKRDTK5CQPQ5SSGEEKG5PTTL5RGGEASERKR
PDSGC5TSKDTKYQ5SVYV5IEEK5DEC
(368-717aa encoded by NM_005618)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.