

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant VP7 protein



Basic Information

Catalog Number:

Ag17616

Size:

50 µg

Form:

Available lyophilized

Species:

Expression Source:

e coli-derived, PET28a, with N-terminal 6[°]His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MDTIAARALTVMRACATLQEARIVLEANVMEILGIAIN
RYNGTLRGV TMRPTSLAQRNEMFFMCLDMMLSAAG
INVGPISPDYTQH MATIGVLATPEIPFTTEAANEIARV
TGETSTWGP ARQPYGFFLETEETFQ PGRWF MRAAQ
AVTAVVCGPDMIQVSLNAGARGDVQQIFQGRNDP
MMIYLVWRR IENFAMAQGNSQQTQAGVTVSVGGV
DMRAGR IIAWDGQAALHVNPTQQNAMVQIQVVF
YISMDKTLNQYPALTAEIFNVYSFRDHTWHGLRTAILN
RTTLPNMLPPIFPNDRDSIL TLLLLSTLADVYTVLRPEF
AIHGVNPM PGLTRAIARAAYV
(1-349 aa encoded by NC_006022)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.