

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human RCC1 protein



Basic Information

Catalog Number:

Ag17826

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MSPKRIAKRRSPADAIPKSKVKVSHRSHSTEPGLVLT
TLGQGDVGQLGLGENVMERKPPALVSIPEDVVQAEA
GGMHTVCLSKSGQVYSFGCNDEGALGRDTSVEGSE
MVPGKVELQEKVVQVSAGDSHTAALDDGRVFLWG
SFRDNNGVI GLLPEMKKSMVPVQVQLDVPVVKVASG
NDHLVMLTADGDLYTLGCGEQQLGRVPELFANRGG
RQGLERLLVPKCVMLKSRGSRGHVRFQDAFCGAYFTF
AISHGHHVYGFGLSNYHQLGTPGTESCIFPQNLTSFK
NSTKSWWGFSGGQHHTVCM DSEGKAYSLGRAEYGRLL
GLGEGAEKSIPTLISRLPAVSSVACGASVGYAVTKDG
RVFAWGMGTNYQLGTGQDEDAWSPVEMMGKQLEN
RVVLSVSSGGQHTVLLVKDKEQS
(1-421aa encoded by BC010067)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.