

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human SP1 protein



## Basic Information

**Catalog Number:**

Ag17828

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

MSDQDHSMDEMTAVVKIEKGVGGNNGGNGGGGA  
FSQARSSSTGSSSTGGGGQESQPSPLALLAATCSRI  
ESPNENSNSQGPSQSGGTGELDLTATQLSQGANG  
WQIISSSSGATPTSKEQSGSSTNGSNGSESSKNRTVS  
GGQYVVAAAPNLQVLTGLPGVMPNIQYQVIPQ  
FQTVDDGQLFAATGAQVQDGSQGIIPGANQ  
QIITNRGSGGNIIAAMPNLLQAVPLQGLANNVLSGQ  
TQYVTNVPVALNGNITLLPVNSVSAATLTPSSQAVTI  
SSSGSQESGSQPVTSGTTISSASLVSSQASSSSFFTN  
ANSYSTTTTTSNMGMNFTTSGSSGTSNGQTPQRV  
SGLQGS DALNIQQNQTSGGSLQAGQQKEGEQNQQ  
TQQQIILIQPLVQGGQALQALQAAPLSGQTFTTQ  
AISQETLQNLQAVPNSGPIIIRTPVGPNGQVSWQ  
TLQLNLQVQNPQAQTITLAPMQGVSLGQTSSSNTT  
LTPIASAASIPAGTVTVNAAQLSSMPGLQITINLSALGT  
SGIQVHPIQGLPLAIANAPGDHGAQLGLHGAGGDGI  
HDDTAGGEEGENSPDAQPQAGRRTREACTCPYCK  
DSEGRGSGDPGKKKQHIHQGGCKVYKTSHLRAHL  
RWHTGERPFMCTWSYCKRFRSDELQRHKRHTGE  
KKFACPECPKRFMRSDHLSKHIKTHQNKGGPGVALS  
VGTLPDLSGAGSESGTATPSALITTNMVAMEAICPE  
GIARLANSGINVMQVADLQISINISGNGF  
(1-785aa encoded by BC062539)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.