

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human FBX05 protein



Basic Information

Catalog Number:

Ag17870

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MSRRPCSCALRPPRCSCSASP AVTAAGRPRPSDSCK
EESSTLSVKMKDFN CNHVSGLKLVK PDDIGRLVSY
TPAYLEGCKDCIKDYERLSCIGSPIVSPRIVELETESKR
LHNKENQHVQQT LNSTNEIEALET SRLYEDSGYSSFSL
QSGLSEHEEGSLEENFGDSLQSC LLQIQSPDQYPNK
NLLPVLHFEKVV CSTLKKNAKRNPKV DREMLKEI IARG
NFR LQNIIGRKMGLECVDI LSELFRRGLRHVLATILAQL
SDMDLINVSKVSTTWKKI LEDDKGAFQLYSKAIQRVT
ENNNKFSPHASTREYVMFRTPLASVQKSA AQTSLKKD
AQTKLSNQGDQKGSTYSRHNEFSEVAKTLKKNESLKA
CIRCNSPAKYDCYLQRATCKREGCGFDYCTKCLCNH
TTKDCSDGKLLKASCKIGPLPGTKKSKK NLRRL
(1-447 aa encoded by BC018905)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.