

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RASAL2 Fusion Protein



Basic Information

Catalog Number:

Ag17975

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

LRFLCPAIMSPSLFNLMQEYPPDDRTSRTLTIKAVIQNL
ANFAKFGNKEEYMAFMNDFLEHEWGGMKRFLLEISN
PDTISNTPGFDGYIDLGRELSVLHSLLEWVVSQLDKAT
VAKLGPLPRVLADITKSLTNPTPIQQQLRRFTEHNSSP
NVSGSLSSGLQKIFEDPTSDLHKLKSPSQDNTDSYFR
GKTLVVQQASSQSMTYSEKDERESSLPNGRSVSLMD
LQDTHAAQVEHASVMLDVPIRLTGSQLSITQVASIKQ
LRETQSTPQSAPQVRRRPLHPALNQPGLQPLSFQNP
VYHLNPIPAMPKASIDSSLENLSTASSRSQNSSEDFK
LSGSPSSSMEDFT

(609-958 aa encoded by BC110611)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.