

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human RAD54L2 protein



Basic Information

Catalog Number:

Ag17995

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

LGGHCLNSSLLVTGQPCGDRHPVLDLRGHKRKLA
TPPAAQESSRRRSRKGHLPAPVQPYEHGYPSVSGGFAMPP
VSLNHNLTTPFTSQAGENSLFMGSTPSYYQLSNLLAD
ARLVFPVTTDPLVPAGPVSSSSTATSVTASNPSFMLN
PSVPGILPSYSLPFSQPLLSEPRMFAPFPSPVLPNSLSR
GMSIYPGYMSPHAGYPAGGLLRSQVPPFDSHEVAEV
GFSSNDEDEKDDDVIEVTGK
(1118-1467 aa encoded by BC024298)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.