

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human C8B protein



## Basic Information

Catalog Number:

Ag18137

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

RNVTEKMASKSGFSFGFKIPGIFELGISSQSDRGKHYIR  
RTKRFSHTKSVFLHARSDLEVAHYKLPKPRSLMLHYEFL  
QRVKRPLEYSYGEYRDLFRDFGTHYITEAVLGGIYEY  
LVMNKEAMERGDYTLNNVHACAKNDFKIGGAIIEEYV  
VSLGVSVGKCRGILNEIKDRNKRD TMVEDLVVLRGG  
ASEHITTLAYQELPTADLMQEWGDAVQYNPAIIKVKV  
EPLYELVTATDFAYSSTVRQNMKQALEEFQKEVSSCH  
CAPCQGNQVPVLKGSRCDCICPVGSQGLACEVSYRK  
NTPIDGKWNWWSNWSGSRRTQRQCNPPPPQN  
GGSPCSGPASETLDCS

(242-591 aa encoded by BC130575)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.