

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

SRGAP2 Fusion Protein



Basic Information

Catalog Number:

Ag18278

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

VQDTEGTVVERSSPKSEI E V I S E P P E E K V T A R A G A S C
P S G G H V A D I Y L A N I N K Q R K R P E S G S I R K T F R S D S H G L S
S S L T D S S S P G V G A S C R P S S Q P I M S Q S L P K E G P D K C S I S
G H G S L N S I S R H S S L K N R L D S P Q I R K T A T A G R S K S F N N H
R P M D P E V I A Q

(784-944 aa encoded by BC132874)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.