

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

IL-27 Fusion Protein



Basic Information

Catalog Number:

Ag18307

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

LLLVQAGVWGFP RPPG R PQLSLQELRREFTVSLHLARK
LLAEV R GQAHRFAESHLP GVNLYLLPLGEQLPDVSLTF
QAWRR L S D P E R L C F I S T T L Q P F H A L L G G L G T Q G R W T N
M E R M Q L W A M R L D L R D L Q R H L R F Q V L A A G F N L P E E E E
E E E E E E E E R K G L L P G A L G S A L Q G P A Q V S W P Q L L S T Y
R L L H S L E L V L S R A V R E L L L S K A G H S V W P L G F P T L S P Q P
(19-243 aa encoded by BC062422)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.