

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human G6PD protein



Basic Information

Catalog Number:
Ag1832

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MAEQVALSRTQVC GILREELFQGD AFHQSDTHIF IIM
GASGDLAKKKIYPTIWWLFRDGLLPENTFIVGYARSRL
TVADIRKQSEPFKATPEEKLKLEDFFARNSYVAGQY
DDAASYQRLNSHMNALHLSQANRLFYALPPTVYEA
VTKNIHESCMSQIGWNRIIVEKPFGRDLQSSDRLSNHI
SSLFREDQIYRIDHYLGKEMVQNLMLVLRANRIFGPIW
NRDNIA CVILTFKEPFGTEGRGGYFDEFGIIRDVMQN
HLLQMLCLVAMEKPASTNSDDVRDEKVKVLCISEVQ
ANNVVLGQYVGNPDGEGEATKGYLDDPTVPRGSTTA
TFAAVVLYVENERWDGVPFILRC GKALNERKAEVRLQ
FHDVAGDIFHQQCKRNELVIRVQPNEAVYTKMMTKK
PGMFFNPEESELDLTYGNRYKNVKLPDAYERLLDVF
GSMHFVRSDELREAWRIFTPLLHQIELEKPKPIPIYIG
SRGPTEADELMKRVGFQYEGTYKWWNPHKL
(1-515aa encoded by BC000337)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.