

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ZNF326 protein



Basic Information

Catalog Number:

Ag18422

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6^{His}.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MGDFGSIHRPGIVVDYQNKSTNVTVAAARGIKRKMM
QPFNKPSGTFIKKPKLAKPMEKISLSKSPTKTDPKNEEE
EKRRIEARREKQRRRREKNSKEYGDGYRMAFTCSFCKF
RTFEEKDIELHLESSSHQETLDHIQKQTKFDKVVMEFL
HECMVNFKFKTSIRKQQTNNQTEVVKIEKDVMEGV
TVDDHMMKIVETVHCSACSVYIPALHSSVQQHLKSPD
HIKGGQAYKEQIKRESVLTATSILNPIVKARYERFVKG
ENSFEIQDHSQDQQIEGD

(207-486 aa encoded by BC090866)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.