

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

HEXA Fusion Protein



Basic Information

Catalog Number:

Ag1854

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MTSSRLWFSLLAAAFAGRATALWPPQNFTSDQ
RYVLYPNNFQFYDVSSAAQPGCSVLDEAFQRYRDL
LFGSGSWPRPYLTGKRHTLEKNVLVSVVTPGCNQLP
TLESVENYTLTINDDQCLLLSETVWGALRGLETFSQLV
WKSAAEGTFFINKTEIEDFPRFPHRGLLDTSRHYLPLSSI
LDTLDVMAYNKLNVFHWHLVDDPSFPYESFTFPELMR
KGSYNPVTHIYTAQDVKEVIEYARLRGIRVLAEFDTPG
HTLSWGPPIGLLTPCYSGSEPSGTFGPVNPSLNNTYE
FMSTFFLEVSSVFPDFYLHLGGDEVDFTCWKSNEIQ
DFMRKKGFGEDFKQLESFYIQLLDIVSSYGKGYVVW
QEVFDNKVKIQPDTIIQVWREDIPVNYMKELELVTKA
GFRALLSAPWYLNRSYGPDWKDFYVVEPLAFEGTPE
QKALVIGGEACMWGEYVDNTNLVPRLWPRAGVAE
RLWSNKLTSDLTFAYERLSHFRCCELLRRGVAQPLNV
GFCEQEFEQT

(1-529 aa encoded by BC018927)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.