

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

C15orf53 Fusion Protein



Basic Information

Catalog Number:

Ag18544

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MEIQGAQEDLGISLSSPRRNHETRPGSKAKGRSSICLQ
ASVWMAAGGKLRRLRASEHLTQGHQQELRDWNLGEDA
SLLFSKSPFGAGKLIQAPAHVFRQCWWQGNAWISCI
TKFDSKRSPEVASSPSYLTVPRRSPLPVFLRPSDRVCV
GGCYLGKSTRRRACQSLLSDPLGVTFPTQTRP
(1-179aa encoded by BC119016)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.