

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human DNASE1L3 protein



Basic Information

Catalog Number:

Ag1873

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MSRELAPLLLLLSIHSALAMRICSFNVRSGESKQEDK
NAMDVIVKVIKRCDIILVMEIKDSNNRICPILMEKLN
SRRGITYNVYVSSRLGRNTYKEQYAFLYKEKLVSVKRSY
HYHDYQDGDADVFSREPFVWFQSPHTAVKDFVIIP
LHTTPETSVKEIDELVEVYTDVKHRWKAENFIFMGDF
NAGCSYVPKKAWKNI RLTPRFVWLI GDQEDTTVKK
STNCAYDRIVLRGQEVSSVVPKSNVDFQKAYKLT
EEEALDVSDHFPVEFKLQSSRAFTNSKKSVTLRKKT
KRS

(1-305 aa encoded by BC015831)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.