

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human TTC26 protein



Basic Information

Catalog Number:

Ag18793

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

CYYKLDYYDVSQEVLAVYLQQIPDSTIALNLKACNHF
RLYNGRAAEAEELKSLMDNASSSFEFAKELIRHNLVVFR
GGEGALQVLPLVDVIPEARLNLVIYLRQDDVQEAY
NLIKDLEPTTPQEYILKGVVNAALGQEMGSRDHMKIA
QQFFQLVGGSAECDTIPGRQCMASCFLLKQFDDV
LIYLSFKSYFYNDIFNFNYAQAKAATGNTSEGEAF
LLIQSEMKNDYIYLSWLARCYIMNKKPRLAWELYLKM
ETSGESFSLQLIANDCYKMGQFYYSAKAFDVLRLDP
NPEYWEGKRGACVGIFQMIAGREPKETLREVLHLLRS
TGNTQVEYMIRIMK

(194-544aa encoded by BC126331)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.