

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

TDRD5 Fusion Protein



Basic Information

Catalog Number:

Ag18896

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSEQERIQECLRKEIRSLLISTKDGLSPQELEKEYLLMV
GNHLPLRILGYRSTMELVLDMPDVVRVCPGAGGTVILK
AIPDESTKGIASLVAKQRSSHKLNRNSMHKGRPSIYSGP
RSHRRVPYRGSVAPILPAVVKSELKDLLALSPVLLSDFE
KAF AKRFGRSFQYMQYGF LSMFEVLNAASDVISVEQ
TRAGSLMLKIKSVTEEKPRGCCPAGKIFTQPFRMKQGS
YSTGFPVAKPCFSQPTSNMEPPKQIMSM EKTSKLN VV
ETSR LNHTEKLNQLENTFKSVIAQIGPGGTISS ELKHKI
KFVVS KFP EGLFISKLLGEYEVIFKEQLSPK KLGFLNVTE
LVGALSD

(1-350 aa encoded by BC130532)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.