

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human UHRF1BP1L protein



Basic Information

Catalog Number:

Ag19085

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

HAHEQDTKMHEIYKGNITPQLNKNTLKTSAAATDVWA
VYFSQFWIDYEGMKSGKGRPI SFVDSFPLSIWICQPT
RYAESQKEPQTCNQVSLNTSQSESSDLA^{GR}LK^{RR}KLL
KEYYSTESEPLTNGGQK^{SS}SDT^{FF}RF^{SP}SS^{SE}ADI^{HLL}
VHVHKHVSMQINHYQY^{LLLL}FL^{HES}L^{LL}SEN^{LR}KDVE
AVTGSPASQTSICIGILL^{RS}AEL^{ALLL}HPVDQANTLKSP
VSESVSPVVPDYLPTENGDFL^{SS}K^{RK}QISR^{DN}R^{IR}SVT
VNHMSDN^{RS}MSV^{DLS}HI^{PL}KD^{PL}LF^{KS}AS^{DT}N^{LQ}KGIS
FMDY^{LS}DKHLGKISE^{DESS}GLVYKSGSGEIG^{SET}SD^{KK}
DSFIQIQVVS

(1089-1440aa encoded by BC127017)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.