

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# USH1C Fusion Protein



## Basic Information

Catalog Number:

Ag1919

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MDRKVAREFRHKVDFLIENDA EKDYLDVLRMYHQT  
MDVAVLVGDLKLVINEPSRLPLFDAIRPLIPLKHQVEYD  
QLTPRRSRKLEVRDLRLHPEGLGSLVRRGGLEFGCGLFI  
SHLIKGGQADSVGLQVGDEIVRINGYSISSCTHEEVIN  
LIRTKKTVSIKVRHIGLIPVKSSPDEPLTWQYVDQFVSE  
SGGVRGSLGSPGNRENKVKVFISLVGSRGLGCSISSG  
PIQKPGIFISHVKPGSLSAEVGLIGDQIVEVNGVDFS  
NLDHKEGRELFMTDRERLAEARQRELQRQELLMQKRL  
AMESNKILQEQQEMERQRKEIAQKAAEENERYRKE  
MEQIVEEEEKFKKQWEEDWGSKEQLLPKTITAEVHP  
VPLRKP KYDQGVPE

(1-392 aa encoded by BC016057)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.