

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

GUCY2F Fusion Protein



Basic Information

Catalog Number:

Ag19229

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

LPEVAARLAIERINRDPSFDLSYSFEYVILNEDCQTSRA
LSSFISHHQMASGFIGPTNPGYCEAASLLGNSWDKGI
FSWACVNYELDNKISYPTFSRTLPSPIRVLVTVMKYFQ
WAHAGVISSDEDIWVHTANRVASALRSHGLPVGVVL
TTGQDSQSMRKALQRIHQADRIIIMCMHSALIGGE
TQMHLECAHDLKMTDGTYYVFPYDALLYSLPYKHTP
YRVLRRNPKLREAYDAVLITVSEQEKFYQAFTEAA
ARGEIPEKLEFDQVSPLFGTIYNSIYFIAQAMNNAMKE
NGQAGAASLVQHSRNMQFHGFNQLMRTDS
(71-398aa encoded by BC156674)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.