

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

A2M Fusion Protein



Basic Information

Catalog Number:

Ag19245

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

KDNGCFRSSGSLNNAIKGGVEDEVTL SAYITIALLEIP
LTVTHPVVRNALFCLES AWKTAQEGDHGSHVYTKAL
LAYAFALAGNQDKRKEVLKSLNEEA VKKDNSVHWER
PQKPKAPVGHFYEPQAPSAEVEEMTSYVLLAYLTAQP
APTSEDLT SATNIVKWITKQNAQGGFSSTQD TVVA
LHALSKYGAATFTRTGKAAQVTIQSSGTFSSKFQVDN
NNRLLQQVSLPELPGEYSMKVTEGECVYLQTS LKYNI
LPEKEEFPFALGVQTLPTCDEPKAHTSFQISLSVSYT
GSRASANMAIVDKMVS GFIPKPTVKMLERSNHVSR
TEVSSNHVLIYLDKVS NQTL SLFFTVLQDVPVRDLKPA
IVKVYDYETDEFAIAEYNAPCSKDLGNA
(1075-1474 aa encoded by BC040071)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.