

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# HIVEP1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag19285

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e. coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

TLQPTGLPSPHHLFSLPLHSQQQSRTPYNNMVPVG  
GIHVVPAGLTYSTFVPLQAGPVQLTIPAVSVVHRTL  
THRNTVTEVSGTTNPA GVAELSSVVPICPIGQIRV  
PGL QNLSTPGLQSLPSLSMETVNIVGLANTNMAPQV  
HPP GLALNAVGLQVLTANPSSQSSPAPQAHIPGLQ  
ILNIAL PTLIPSVSQVAVDAQGAPEMPASQSIKACET  
QPKQTS VASANQVSRTESPQGLPTVQRENAK  
KVLNPPAPAGD HARLDGLSKMDTEKAASAN  
HVKPKPELTSIQGQPAST SQPLLKAHSEVFT  
KPSGQQTLSPDRQVPRPTALPRRQ PTVHFSD  
VSSDDEDLVIAT (2366-2718aa encoded by BC140816)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.