

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# FLRT2 Fusion Protein



## Basic Information

Catalog Number:

Ag19435

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6<sup>His</sup>.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

PSTASPTTQPPTLSIPNPSRSYTPPTPTTSKLPDWD  
GRERVTPPISERIQLSIHFVNDTSIQVSWLSLFTVMAY  
KLTWVKMGHSLVGGIVQERIVSGEKQHLSLVNLEPRS  
TYRICLVPLDAFNRAVEDTICSEATTHASYLNNGSNT  
ASSHEQTT

(374-533 aa encoded by BC126249)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.