

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human CD41/Integrin alpha 2b protein



## Basic Information

### Catalog Number:

Ag19459

### Size:

50 µg

### Form:

Available lyophilized

### Species:

human

### Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

### Biological Activity:

Not tested

### Endotoxin Level:

Please contact the lab for more information

### Peptide Sequence:

AELAVHLPQGAHYMRALSNVEGFERLICNQKKENETR  
VVLCELGNPMKKNAQIGIAMLVSVGNLEEAGESVSF  
QLQIRSKNSQNPNSKIVLLDVPVRAEAQVELRGNSFP  
ASLVVAAEEGEREQNSLDSWGPKEHTYELHNNGPG  
TVNGLHLSIHLPGQSQPSDLLYLDIQPQGGLQCFQP  
PVNPLKVDWGLPIPSPIHPAHHKRDRRQIFL  
(678-894 aa encoded by BC126442)

## Selected Validation Data