

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

BCMO1 Fusion Protein



Basic Information

Catalog Number:

Ag19519

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

FKIPATVPEGKKGKSPWKHTEVFCISPSRLLSPSY
HSFGVTENYVIFLEQPFRLDILKMATA YIRMSWASCL
AFHREEKTYIHIDQRTRQPVQTKFYTDAMVVFHHVN
AYEEDGCIVFDVIA YEDNSLYQLFYLANLNQDFKENS
RLTSVPTLRRFAVPLHVDKNAEVGTNLIKVASTTATAL
KEEDGQVYQCPEFLYEGLELPRVNYAHNGKQYRYVF
ATGVQWSP IPTKIIKYDILT KSSLKWREDDCWPAEPLF
VPAPGAKDEDDGVILSAIVSTDPQKLPFLIILDAKSFTE
LARASVDVDMHMDLHGLFITDMDWDTKKQAASEEQ
RDRASDCHGAPLT

(199-547aa encoded by BC126210)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.