

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human PLEKHG7 protein



Basic Information

Catalog Number:

Ag19560

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MIFMNTLRYLQTHEYLLDVLWRLFANLEELTQTSLGF
VNSLFGIHKDYVDASEISSSLDFISVLTKYFRGSLCQSH
QTYCLNYSAAIFYLES LRQRDDFGIY LKWCEQNEQCR
RLHVPELLVAPLQRLTRYPLLLKNIWKRSMDSAEKIMY
SIKEKVEKSI RDLEGGV KWLDFNQKFRYLQEIIVWPPL
WDRDKRFFIPECLKHIFKEHMAENISPTSRLHLYEGKL
TLAESTRFLDVYLFNFDFLLVTKCKNKKLGGSDPGL
MCPSLTPELQAVIKEGGCTVLDQPIPLDRLVVKSIEPL
HVSVFLGRNAFLIQHENRYRQCIAAFLLQAQTENIKK
TWMAQITTAISCFTKSQETKKISLFTLPAESSEI
(1-379aa encoded by BC137086)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.