

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human STK32C protein



Basic Information

Catalog Number:

Ag1960

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

KKRMGSSMSAATARRPVFDDKEDVNFDFHFIILRAIGK
GSFGKVCIVQKRDEKMYAMKYMKNQQCIEERDEVNRN
VFRELEILQIEHVFLVNLWYSFQDEEDMFVVDLL
GGDLRYHLQQNVQFSEDTVRLYICEMALDYLRGQ
HIIHRDVKPDNILLDERGHAHLTDFNIATIIKDGERTATA
LAGTKPYMAPEIFHSFVNGGTGYSFEVDWWSVGVGM
AYELLRGWRPYDIHSSNAVESLVQLFSTVSVQYVPTW
SKEMVALLRKLLTVNPEHRLSSLQDVQAAPALAGVLW
DHLSEKRVEPGFVPNKGRLHCDPTFELEEMILESRLPH
KKKKRLAKNKS RDNSRDSSQSENDYLQDCLDAIQQDF
VIFNREKLKRSQDLPREPLPAPESRDAAEPVEDEAERS
ALPMC GPICPSAGSG

(1-422 aa encoded by BC015792)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.