

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human LAMC3- Specific protein



## Basic Information

Catalog Number:

Ag19613

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

SLSEGGCRPCTCNPAGSLDTC DPRSGRPC KENVEG  
NLCDRCRPGTFNLQPHNPAGC S S C F C Y G H S K V C A S T  
A Q F Q V H H I L S D F H Q G A E G W W A R S V G G S E H P P Q W S  
P N G V L L S P E D E E E L T A P E K F L G D Q R F S Y G Q P L I L T F R V P  
P G D S P L P V Q L R L E G T G L A L S L R H S S L S G P Q D A G H P R E V  
E L R F H L Q E T S E D V A P P L P P F H F Q R L L A N L T S L R L R V S P G  
P S P A G P V F L T E V R L T S A R P G L S P P A S W V E I C S C P T G Y T  
G Q F C E S C A P G Y K R E M P Q G G P Y A S C V P C T C N Q H G T C  
D P N T G I C V C S H H T E G P S C E R C L P G F Y G N P F A G Q A D D  
C Q

(421-753 aa encoded by NM\_006059)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.