

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ZNF639 protein



Basic Information

Catalog Number:

Ag19854

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

EDVPIAVEVHAISEDYDIETENSSSLSQDQTDDEPP
AKLCKILDKSQALNVTAQQKWPLLRANSSGLYKCELC
EFNSKYFSDLKQHMLKHKRTDSNVCRVKESFSTNM
LLIEHAKLHEEDPYICKYCDYKTVIFENLSQHIADTHFS
DHLYWCEQCQDVQFSSSELYLHFQEHSCDEQYLCQF
CEHETNDPEDLHSHVVNEHACKLIELSDKYNNGEHG
QYSLLSKITFDCKNFFVCQVCGRSRLHTNVNRHVA
IEHTKIFPHVCDDCGKGFSSMLEYCKHLNSHLEGIYL
CQYCEYSTGQIEDLKIHLDFKHSADLPHKCSDCMLRF
GNERELISHLPVHETT

(136-485 aa encoded by BC026181)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.