

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human CASP8AP2 protein



## Basic Information

Catalog Number:

Ag19912

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

KDASDKMGHSDEVADECFLHQVWETKVPESIEELP  
SMEEISHSVGEHLPNTYVDLTQDPVTETKNLGEFIEVT  
VLHIDQLGCSGGNINQSAQILDNSLQADTVGAFIDLT  
QDASSEAKSEGNPALAVEDLGCGVIQVDEDNCKEEK  
AQVANRPLKCIVEETYIDLTTESSCEVKKDELKSEP  
GSNCDNSELPGTLHNSHKRRNISDLNHPHKQRKET  
DLTNKEKTKKPTQDSCENTEAHQKASKKKAPPVTKD  
PSSLKATPGIKDSSAALATSTSLAKNVIKKKGIIILWT  
RNDREILLECQKRGPSFKTFAYLAAKLDKNPNQVSR  
FQQLMKLFEKSKCR

(1631-1981 aa encoded by BC150274)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.