

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# SMTN Fusion Protein



## Basic Information

**Catalog Number:**

Ag19983

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

KTTFTIEIKDGRGQASTGRVLLPTGNQRAELTLGLRAP  
PTLLSTSSGGKSTITRVNSPGLARLGSVTHVTSFSHA  
PPSSRGGCSIKMEPEPAEPLAAAVEAANGAEQTRVN  
KAPEGRSPLSAEELMTIEDEGVLDKMLDQSTDFEERKL  
IRAAALRELQRKRQDQDKERERRLQEARGRPGEGRGN  
TATETTTRHSQLAADGSAVSTVTKTERLVHSDGTRT  
ARTTTVESSFVRRSENGSGSTMMQTKTFSSSSSSKMM  
GSIFDREDQASPRAGSLAALEKRQAEKKKELMKAQSL  
PKTSASQARKAMIEKLEKEGAAGSPGGPRAAVQRSTS  
FGVNPANSIKQMLLDWCRA  
(458-811 aa encoded by BC034237)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.