

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human FAM134A protein



Basic Information

Catalog Number:

Ag20025

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

DYSMKA E A N A L H H K H D K R K R Q G K N A P P G G D E P L A E T
E S E S E A E L A G F S P V V D V K K T A L A I T D S E L S D E E A S I L
E S G G F S V S R A T T P Q L T D V S E D L D Q Q S L P S E P E E T L S R
D L G E G E E G E L A P P E D L L G R P Q A L S R Q A L D S E E E E E D V
A A K E T L L R L S S P L H F V N T H F N G A G S P P D G V K C S P G G P
V E T L S P E T V S G G L T A L P G T L S P P L C L V G S D P A P S P S I L P
P V P Q D S P Q P L P A P E E E E A L T T E D F E L L D Q G E L E Q L N A
E L G L E P E T P P K P P D A P P L G P D I H S L V Q S D Q E A Q A V A E
P

(244-543 aa encoded by BC002420)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.