FOR IN VITRO RESEARCH USE ONLY. NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human C11orf48 protein



Basic Information

Catalog Number:

Ag20028

Size:

50 μg Form:

Available lyophilized

Species: human

Expression Source:

e coli.-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

LPHEELLLTDGEEEDAEAFFQDQSEEPGWAWSPQD
PRSPLRTFNAGLSWGQDQDEEDACWILEDTACLEAT
NHCPFWDSTGSRVCRSGFVEYSHLLPPNSFEGAEEEA
VQTPAGVESGAASEAPGGRGCDRPRADHAAPPQEA
GVQCTCQHYTVREEAQKTPPADPACPEREDSHGSG

(77-289 aa encoded by BC001434)

Selected Validation Data