

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human KRT4 protein



## Basic Information

Catalog Number:

Ag20210

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MIARQQCVRGGPRGFSCGSAIVGGGKRGAFSSVSMS  
GGAGRCSSGGFGSRSLYNLRGNKSISMSVAGSRQGA  
CFGGAGGGFTGGFGAGGFGAGFTGGFGGGFGGSF  
SGKGGPGFPVCPAGGIQEVNTINQSLLTPLHVEIDPEIQ  
KVRTEEREQIKLLNNKFASFIDKVQFLEQQNKVLETK  
WNLLQQQTTTTSSKNLEPLFETYLSVLRKQLDTLGNDK  
GRLQSELKTMQDSVEDFKTKYEEIINKRTAAENDFVV  
LKKDVAAYLNKVELEAKVDSLNDIENFLKVLDAELS  
QMQTHVSDTSVVLMDNNRNLDSLIIAEVRAQYEEI  
AQRSKAEAEALYQTKVQQLQISVDQHGDNLKNTKSE  
IAELNRMIRLRAEIEINIKKQCQTLQVSVADAEQRGE  
NALKDAHSKRVELEAALQQAKEELARMLREYQELMSV  
KLALDIEIATYRKLLGEEYRMSGECQSAVSISVSVSGST  
STGGISGGLGSGGFLSSGFGSGSGFGFGGSVSG  
SSSSKIISTTTLNKR

(1-534 aa encoded by BC042174)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.