

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human Ubiquitin 2 protein



Basic Information

Catalog Number:

Ag20251

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MAENGESSGPPPSRGPAQAQGSAAAPAEPKIKVTV
KTPKEKEEFVAVPENS SVQFKEAISKRFKSTQDQLVLI
FAGKILKDQDTLIQHGIHDGLTVHLVIKSNRPQGQS
TQPSNAAGTNTTSASTPRSNSTPISTNSNPFGLGSLG
GLAGLSSLGLSSTNFSELQSQMQQLMASPEMMIQI
MENPFVQSMLSNPDLMRQLIMANPQMQLIQRNPEI
SHLLNPDIMRQTLEIARNPAMMQEMMRNQDLALSN
LESIPGGYNALRRMYTDIQEPLNAAQEQFGGNPFAS
VGSSSSSGEGTQPSRTENRDPLNPWAPPATQSSA
TTSTTTSTGSGSGNSSNATGNTVAAANYVASIFSTP
GMQSLQKQITENPQLIQNMLSAPYMRSMQSLSQNP
DLAAQMMNLSPFLTANPQLQEQMRPQLPAFLQMQ
NPDTLSAMSNPRAMQALMQIQQLQLATEAPGLIP
SFTPGVGVGLGTAIGPVGPVTPIGPIGPIVPFTPIGPI
GPIGPTGPAAPPSTGSGGPTGPTVSSAAPSETTSPT
SESGPNQQFIQQMVQALAGANAPQLPNPEVRFQQ
QLEQLNAMGFLNREANLQALATGGDINAAIERLLGS
QPS

(1-624 aa encoded by BC069237)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.