

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human NFkB1 protein



## Basic Information

Catalog Number:

Ag20361

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

LAGCLLLEGDHVDSTTYDGTTPHIAAGRGSTRLAAL  
LKAAAGADPLVENFEPLYLDDSWENAGEDEGVVPGT  
TPLDMATSWQVFDILNGKPYEFTSDDLLAQGDMK  
QLAEDVKLQLYKLEIPDPDKNWATLAQKLGILNNA  
FRLSPAPSKTLMNDNYEVSGGTVRELVEALRQMGYTEA  
IEVIQAASSPVKTTSAHSLPLSPASTRQQIDELRSD  
SVCDSGVETSFRKLSFTESLTSGASLLTLNKMPHDYGO  
EGPLEGKI

(701-970 aa encoded by BC051765)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.