

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human MMAA protein



## Basic Information

**Catalog Number:**

Ag20370

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET30a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

TTLKDHT EGLSDKEQRFVDKLYTGLIQGQRACLA EAIT  
LVESTHSRKKELAQVLLQKVLVLYHREQEQSNKGKPLAF  
RVGLSGPPGAGKSTFIEYFGKMLTERGHKLSVLA VDPS  
SCTSGGSLGDKTRMTELSRDMNAYIRPSPTRGTLGGV  
TRTTNEAILLCEGAGYDII LIETVGVGQSEFAVADMVD  
MFVLLPPAGGDELQGIKRGIIEMADLVAVTKSDGDLI  
VPARRIQAEYVSALKLLRKRQVWPKVIRISARSGEG  
ISEMWDKMKDFQDLMLASGELTAKRRKQKQVWMW  
NLIQESVLEHFRTHPTVREQIPLEQKVLIGALSPGLAA  
DFLLKAFKSRD

(69-418 aa encoded by BC101179)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.