

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human BCL3 protein



Basic Information

Catalog Number:

Ag21040

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET30a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

NLPTPLYPM MCPMEHPLSADIAMATRADEGDGDTPLHI
AVVQGNLPAVHRLVNLFQQGGRELDIYNNLRQTPLHL
AVITTLPSVVRLVLTAGASPMALDRHGQTAHLACEH
RSPTCLRALLDSAAPGTLDEARNYDGLTALHAVNTE
CQETVQLLLERGADIDAVDIKSGRSPLIHAVENNSLSM
VQLLLQHGANVNAQMYSGSSALHSASGRGLPLVRTL
VRSGADSSLNCHNDTPLMVARSRVIDILRGKATRP
ASTSQPDPSDRSANTSPSSSRLSSNGLLSASPSSP
SQSPPRDPPGFMAPPNFFLPSPSPAFLPFAGVLRGP
GRPVPPSPAPGGS

(97-446 aa encoded by BC064993)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.