

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

TLE1 Fusion Protein



Basic Information

Catalog Number:

Ag21327

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

DPPPHMRVPTIPPNLGIPGGKPAYSFHVTADGQMQ
PVPFPDALIGPGIPRHARQINTLNHGEVCAVTISNP
TRHVYTGKGCVKVWDISHPGNKSPVSLDCLNRDN
YIRSKLLPDGCTLIVGGEASTLSIWDLAAPTPIKAEL
TSSAPACYALAI SPDSKVCFCSDGNIAVWDLHNQ
TLVRQFQGHTDGASCIDISNDGTKLWTGGLDNTVRS
WDLREGRQLQHDFTSQIFSLGYCPTGEWLA VGMES
SNVEVLHV NKPDKYQLHLHESCVLSLKFAYCGKWFVS
TGKDNLLNAWRTPYGASIFQSKESSVLSCDISVDDK
YIVTGSGDKKATVYEVY

(422-770 aa encoded by BC010100)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.