

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human ITPKB protein



Basic Information

Catalog Number:

Ag21362

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MAVYCYALNSLVMNSANEMKSGGGPGPSGSETPPP
PRRAVLSPGSVFSPGRGASFLFPPAESLSPPEPRSPGG
WRSGRRLNSSSGSGSGSSVSSPSWAGRLRGDR
QQVVAAGTLPSPGPEEAKRKLRLQRELQNVQVNQK
VGMFEAHIQAQSSAIQAPRSPRLGRAHSPSPCFR55
SQPPGRVLVQGARSEERRTKSWGEQCPE5TSGTDSGR
KGGPSLCSSQVKKGMPPLPGRAAPTGSEAQG5SAFV
RMEKGIPASPRCGSPTAMEIDKRG5PTPGTR5CLAPSL
GLFGASLT5MATEVAARVT5TGPHRPQDLALTEPSGRA
RELEDLQPPEALVERQ5QFL

(1-350 aa encoded by BC015009)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.