

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

GSK3B Fusion Protein



Basic Information

Catalog Number:

Ag21438

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET30a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSGRPRTTSFAESCKPVQQPSAFGSMKVS RDKDGSK
VTTVVATPGQGPD RPQEVSYDTKVI GNGSFGVVYQ
AKLCDSGELVAIKKVLQDKRFKNRELQIMRKL DHCNIV
RLRYFFYSSGEKKDEVYLNVLVDYVPETVYRVARHYSR
AKQTL PVIYVKLYMYQLFRSLAYIHSFGICHRDIKPQNL
LLDPDTAVLKLCDFGSAKQLVRGEPNVSYICSRYYRAP
ELIFGATDYTSSIDVWSAGCVLAELLGQPIFPGD SGV
DQLVEI IKVLGTP TREQIREMNP NYTEFKFPQIK AHPW
TKDSSGTGHFTSGVRVFRPRTPEAIALCSR LLEYTPT
ARLT PLEACAHSFFDEL RDPNVKLPNGRDT PALFNFTT
QELSSNPPLATILIPPHARIQAAA STPTNATAASDANT
GDRGQTNNAASASASNST
(1-433 aa encoded by BC000251)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.