

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant rat Angptl4 protein



Basic Information

Catalog Number:

Ag21514

Size:

50 µg

Form:

Available lyophilized

Species:

rat

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MRCAPTAGAALVLC AATAGLLSAQGRPAQPEPPRFA
SWDEMNL LAHGLLQLGHGLREHVERTRGQLGALERR
MAACGNACQGP KGTDPKDRVPEGQAPETLQSLQTQ
LKAQNSKIQLFQKVAQQQRYLSKQNLRIQNLQSQI
DLLTPTHLDNGVDKTSRGKRLPKMAQLI GLTPNATRLH
RPPRDCQELFQEGERHSGLFQIQPLGSPFFLVNCEMT
SDGGWTVIQRR LNGSVDFNQSW EAYKDGFGDPQGE
FWLGLEKMH SITGDRGSQLAVQLQDWDGNAKLLQFP
IHLGGEDTAYS LQLTEPTANELGATNVSPNGLSLPFST
WDQDHDLRGDLNCAKSLSGGWFFGTCSHSNLNGQ
YFHSIPRQRQQRKKGIFWKTWKGRYYPLQATTLLIQP
MEATAAS

(1-405 aa encoded by BC078944)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.