

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

KIAA0020 Fusion Protein



Basic Information

Catalog Number:

Ag22011

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET30a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

YKSADHPTLDKVLQPEKLEIMDEMKGILTPMAQKE
AVIKHSLVHKVFLDFFTYAPPKLRSEMIEAIREAVVYLA
HTHDGARVAMHCLWHGTPKDRKVIVKTMKTYVEKV
ANGQYSHLVLLAAFDCIDDTKLVKQIIISEISSLPSIVN
DKYGRKVLVLLSPRDPAAHTVREIIEVLQKGDGNAHSK
KDTEVRRRELLESISPALLSYLQEHAEVVLDKSACVL
VSDILGSATGDVQPTMNAIASLAATGLHPGGKDGELH
IAEHPAGHLVLKWLIEQDKKMKENGREGCFAKTLVEH
VGMKNLKSASVNRGAILSSLLQSCDLEVANKVKAA
LKSLIPTLEKTKSTSKGIEILLEKLST
(283-648 aa encoded by BC016137)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.