

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human MEX3D protein



Basic Information

Catalog Number:

Ag22080

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET30a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

VC L D L L G A A A S L W A K T P N Q G R R P P T A T A G L R G D T A L
G A P S A P E A F Y A G S R G G P S V P D P G P A S P Y S G S G N G G F
A F G A E G P G A P V G T A A P D D C D F G F D F D L A L D L T V P A
A A T I W A P F E R A A P L P A F S G C S T V N G A P G P P A A G A R R
S S G A G T P R H S P T L P E P G G L R L E L P L S R R G A P D P V G A L S
W R P P Q G P V S F P G G A A F S T A T S L P S S P A A A A C A P L D S
G A S E N S R K P P S A S S A P A L A R E C V V C A E G E V M
(184-432 aa encoded by BC113384)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.