

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

WDR16 Fusion Protein



Basic Information

Catalog Number:

Ag22166

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET30a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

GEGHQFLVGT EESH IYRV SFTDFK ETLIATCHFD AVK
DIVFPFGTAE LFATCAKKDIRVWHTSSNRELLRITV PN
MTCHGIDFMRD GKSII SAWNDGKIRAFAPETGR LMYV
INNAHRIGVTAI ATTS DCKRVISGGGEGEVRVWQIGC
QTQKLEEA LKEHKSSV SCIRVKRNNEECVTASD GTC
IIWDLVRLRRNQ MILANTLFQCVCYHPPEFQIITSGTD
RKIAYWEVFDGT VIRELEGLSGSINGMDITQEGVHF
VTGGNDHLVKVWDYNEGEVTHVGVGHSGNITRIRIS
PGNQYIVSVSADGAILRWKYPYTS

(300 - 620 aa encoded by BC025392)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.