

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# TPRX1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag23015

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PET30a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

PGCRAQKGI PAALSPGPGPI PAPI PGP AQIPG PLPGSIP  
GPI PGPAQIPSP I PAPI PGPI SGVQIPG PFRGPI PGPI S  
GPAPI PGPI SGPFSGPNPGPI PGPNPGPI PGPI SGPI PG  
PISVPI PGPI PGPI SGPI SGPNPGPI PGPI SGPNPG  
PI PGPI SGPNPGLI PGPI PGPI SGPGPI I GPI SPAQIPGP  
GRLQGGPI LSPGRMRSPGSLPGLAPILGPGSGPGSGS  
VPAPI PGPSLPAPAPLWPQSPDASDFLPDTQLFPHF  
TELLPLDPLEGS SVTMTSQYQEGDDSMGKKHSGSQ  
PQEEGGSVNENHSGPRLLDL  
(80-411 aa encoded by BC137501)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.