

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

WAS Fusion Protein



Basic Information

Catalog Number:

Ag23331

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e. coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MSGGPMGGRPGGRGAPAVQQNIPSTLLQDHENQRL
FEMLGRKCLTLATAVVQLYLALPPGAEHWTKEHCGA
VCFVKDNPQKSYFIRLYGLQAGRLLWEQELYSQLVYS
TPTPFHTFAGDDCQAGLNFADEDEAQAFRALVQEK
IQKRNQRQSGDRRQLPPPPTPANEERRGGLPPLHP
GGDQGGPPVGPLSLGLATVDIQNPDISSRYRGLPAP
GPSPADKKRSGKKKISKADIGAPSGFKHVSHVSWDP
QNGFDVNNLDPDLRSLFSRAGISEAQLTDAETSKLIYD
FIEDQGGLEAVRQEMRRQEPLPPPPPSRGGNQLPRP
PIVGGNKGRSGPLPPVPLGIAPPPPTPRGPPPPGRGGP
PPPPPPATGRSGPLPPPPPGAGGPPMPPPPPP
(1-400 aa encoded by BC012738)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.